Test Solutions for Medical Device Manufacturers



Mastering Verification and Validation for Multi-Vital Sign Testing in R&D

Speaker: Stephan Syue

Outline

Time	Contents
2.5 min.	 Vital Sign Testing Concepts The Correlation between Vital Signs – Taking Patient Monitors as An Example Multiparameter Simulation Demonstration – ECG, PPG, Impedance Respiration, and NIBP Signals
23 min.	 ECG & PPG Testing Brief Introduction of ECG Compliance Tests PPG Principle for SpO₂ Measurement Transmittance SpO₂ Module PPG-2TF-660 Demonstration
18 min.	 Respiration & NIBP Testing Optimized Respiration Simulation and Raw Data Playback Brief Introduction of NIBP Compliance Tests NIBP Simulator BPA700 Demonstration
4.5 min.	Automated Testing Solution • Sample Code of ECG and NIBP Compliance Tests • User Testimonial – ECG Automated Testing
10 min.	Conclusions and Q&A



(1) Vital Sign Testing Concepts

Vital Sign Correlation Powered By Multiparameter Simulation

Next-generation Medical System



More Accurate and Effortless Diagnosis : Vitals Correlation+Al

Data-says-it-all Medical Future





(2) ECG & PPG Testing





(2-1) ECG Testing SECG 5.0 AIO Multi Vital Sign Simulator





- 2-1-2 Test Method
- ► 2-1-3 Test Example

Six Essentials of ECG Standard Test Items

2-1-1 IEC 60601-2-27 Test Requirements

For example: IEC 60601-2-27:2011 Patient Monitoring ECG

201.7.9.2.9.101 Additional instructions for use

201.7.9.2.9.101 b) 4) Heart rate meter accuracy and response to irregular rhythm 201.7.9.2.9.101 b) 5) Response time of heart rate meter to change in heart rate 201.7.9.2.9.101 B) 6) Time to alarm for tachycardia

201.12.1.101 ESSENTIAL PERFORMANCE of ME EQUIPMENT

201.12.1.101.1 Accuracy of signal reproduction 201.12.1.101.2 Input dynamic range and differential offset voltage 201.12.1.101.3 Input impedance 201.12.1.101.4 Input NOISE 201.12.1.101.5 Multichannel crosstalk 201.12.1.101.6 GAIN control and stability 201.12.1.101.7 Sweep speed 201.12.1.101.8 Frequency and impulse response 201.12.1.101.9 GAIN INDICATOR 201.12.1.101.10 Common mode rejection 201.12.1.101.11 Baseline reset 201.12.1.101.12 Pacemaker pulse display capability 201.12.1.101.13 Rejection of pacemaker pulses 201.12.1.101.14 Synchronizing pulse for cardioversion 201.12.1.101.15 Heart rate range, accuracy, and QRS detection range 201.12.1.101.16 Channel height and aspect ratio 201.12.1.101.17 Tall T-wave rejection capability

ALGORITHM PERFORMANCE

Verify with PQRST Complex, Clinical Database, etc.

ELECTRICAL PERFORMANCE

- Amplitude Accuracy
- Dynamic Range
- Input Impedance
- Frequency Response

Noise

Verify with Basic Waveforms

Advanced Algorithm Testing

Machine Learning as Looking up to Annotation



ECG Test Capabilities powered by Hassle-free HW and GUI



2-1-2 Test Method

Test Example

201.12.1.101.2 Input dynamic range and differential offset voltage

- Purpose
 - With a d.c. offset voltage in the range of $\pm 300 \text{ mV}$ and differential input signal voltages of $\pm 5 \text{ mV}$ that vary at rates up to 320 mV/s, when applied to any LEAD WIRE, the timevarying output signal amplitude shall not change by more than $\pm 10 \%$ over the specified range of d.c. offset.
- METHOD
 - Open S1, close S and S2 and S4 at position B (no DC).
 - Apply a 16 Hz triangular or sinusoidal signal to any LEAD
 WIRE with all other LEAD WIRES connected to the N (RL)
 LEAD WIRE (P2)as defined in Table 201.103.
 - Set the GAIN to 10 mm/mV and sweep speed to 25 mm/s.
 - Adjust the signal generator so that the applied input signal produces an output amplitude of 80% of the full scale channel height. Record the amplitude of this output signal.
 - Set S4 to position A (with DC)to apply a d.c. offset voltage of +300 mV.
 - Repeat this test for a d.c. offset voltage of -300 mV by changing the position of S3.
 - Repeat the test for each LEAD WIRE until all combinations of LEAD WIRES have been tested as defined in Table 201.103.
 - Repeat the test for each PERMANENT DISPLAY and NON-PERMANENT DISPLAY.
- RESULTS
 - Verify that the output signal amplitude is within ±10% of the previously recorded amplitude

(Triangle wave: SR=4fVp=4x 16 x 5 mV=320 mV/s)



EUT Display & Measurement :



Test Example

201.12.1.101.2 Input dynamic range and differential offset voltage



Advanced Algorithm Testing

Machine Learning as Looking up to Annotation





(2-2) PPG Testing PPG-2TF-660 Transmittance SpO₂ Module



- 2-2-1 PPG Technology for SpO₂ Measurement
- 2-2-2 Test Method
- ► 2-2-3 Operational Principle
- 2-2-4 Test Method

Three Essentials of PPG for SpO₂



2-2-1 PPG Technology for SpO₂ Measurement



© 2025 WhaleTeq All Rights Reserve

Verify Readings with Recorded Raw Data or Customized Waveform Replayed by SECG 5.0 AIO

Advanced Algorithm Testing

Machine Learning as Looking up to Annotation



WhaleTeqSECG50		×	Pulse Oximeter Signa	I-in
WhaleTeq Single Channel Signal Type	Load Raw Files	2025/1/8 04:27PM	2	
ВРМ	Channel: PPG Red C:\Users\Stephan\AppData\Local\Programs\WhaleTeq SECG 5.0\User\Raw Data\ Channel [] [] []		Recording Replay	PPG 3
	[mgnml-waveform-database-1.0.0] lightsleep_ppg_200.txt s1_walk_gz_Clinical PPG X2.txt			Measurement
Sampling Load Raw File	Output Setting AC 12.50 mV DC 625 mV 0 + 0.121 x ([data] + 0.000)	Ambient Light OFF	Record wit Annota	ding h ation
Autoritation and a constraint of the second	OK Cancel	9.0 10.0 STOP		

Operational Principle

Synchronize DUT Sampling Frequency



2-2-3 Operational Principle

PPG-2TF-660 Functions

R&D-purposed Transmittance PPG Simulator



Load User-defined Waveform to Replay for Algorithm Enhancement



Derive DUT-tailored Reference R Curve to Check Repeatability and Reproducibility



SpO ₂ (%)	SN 1 Golden	SN 2	SN 3	SN 4
95	<mark>95</mark> 95	94	95	96
90	89			
85	85			
80	81			

2 - 2 - 4

Test Method

Add-on Signals: Ambient Noises, Finger type-specific Transmission Rates & Respiration







- Flexible ECG, PPG, Impedance RESP and NIBP Simulations At Once enhances analysis of vital sign correlation by patient monitors.
- Standard Assistant Software breaks down and streamlines ECG medical standard tests into Pass Criteria, Preparation, Test Parameter per Condition and Auto-verdict.
- Independent R/IR-AC and R/IR-DC Settings of PPG Simulation portrays various patient conditions.
- DUT-tailored R Curve helps benchmark new product with the existing one and also assures repeatability and reproducibility.
- Replaying clinical ECG and / or PPG database the analogue way cost-effectively enhances algorithm testing.



(3) Respiration & NIBP Testing





(3-1) Respiration Testing

Respiration Module



- ▶ 3-1-1 Respiration Module Introduction
- ▶ 3-1-2 Test Method

Advanced Respiration Parameters

Both Impedance and Modulation Methods Available





Advanced Algorithm Testing

Machine Learning as Looking up to Annotation







(3-2) NIBP Testing BPA700

NIBP Simulator & Analyzer



3-2-1 IEC 80601-2-30 Test Requirements

3-2-2 Test Method

- 3-2-3 IEC 80601-2-30 Test Requirements
- ► 3-2-4 Test Example
- ▶ 3-2-5 IEC 80601-2-30 Test Requirements

▶ 3-2-6 Test Example

▶ 3-2-7 BPA700 Specifications

Six Essentials of NIBP Standard Test Items

IEC 80601-2-30

 201.11.8: Interruption of the power supply/supply mains to the equipment. 201.11.8.101 Switching off 201.11.8.102 Interruption of the supply mains
 201.12.1: Accuracy of controls and instruments. 201.12.1.101 Measuring and display ranges 201.12.1.102 Limits of the error of the manometer from environmental conditions 201.12.1.103 Nominal blood pressure indication range 201.12.1.104 Maximum pressure in normal condition 201.12.1.105 Maximum pressure in single fault condition 201.12.1.106 Manometer test mode

201.12.1.107 Reproducibility of the blood pressure determination

3 201.101.2: Pressurization.

4 201.104: Maximum inflating time.

(5) 201.105: Automatic cycling mode.

201.105.1 Long-term automatic mode 201.105.2 Short-term automatic mode 201.105.3 Self-measurement automatic mode 201.105.3.1 General 201.105.3.2 Normal condition 201.105.3.3 Single fault condition 1 Deflation Time

(2) Manometer Accuracy & Sys/Dia Repeatability

③ Protection Mechanism at Max. Pressure

(4) Cuff and Bladder Robustness

(5) Max. Inflating Time

6 Sys/Dia Accuracy

Compliance **#** Competence



Friendly GUI BPA700







201.12.1.106 * Manometer test mode

The AUTOMATED SPHYGMOMANOMETER shall have a manometer test mode that permits static pressure measurement over at least the NOMINAL BLOOD PRESSURE indication range (see 201.12.1.103). This mode shall not be available in NORMAL USE but restricted to SERVICE PERSONNEL.

EXAMPLE 1 A port for connection to a pressure source so that the pressure can be measured by the AUTOMATED SPHYGMOMANOMETER in a test mode.

Manometer of NIBP Monitor

EXAMPLE 2 A port for connection to a reference manometer that can be pressurized by the AUTOMATED SPHYGMOMANOMETER in a test mode. Pressure Source of NIBP Monitor

NOTE This mode can be used to verify manometer pressure accuracy.

The technical description shall include a test method that can be used to verify the calibration of the AUTOMATED SPHYGMOMANOMETER.

Compliance is checked by inspection and functional testing.

Static Pressure Test



against Manometer of NIBP Monitor



Static Pressure Test







BPA700 Standard Assistant



201.12.1.106 Manometer test mode

Assistant			— C	×						
Standard:	IEC 80601-2-30:2018	·		=						
	201.12.1.106 Manometer test mode	· · · · · · · · · · · · · · · · · · ·	·							
Pass Criteria:	Over at least the nominal blood pressure available in normal use, but restricted to s	indication range (20/ i 10 mmHg and 40 service personnel.	//230 mmHg). This mode shall not	pe						
Preparation:	1. Connect the automated sphygmoman 2. Connect air reservoir to port air reservo	ometer to DUT. ir (Other mode: 500ml, Neonatal mode	: 100ml).							
	Note: It is recommended not to drasticall mmHg.	y change blood pressure values, e.g., st	atic pressure from 300 mmHg to 1	0 Pressure G	iraph					
Test Parameter:	Rated range: Min 0 🚔 Max 260	Measure Time 5 second	nd(s)	3 00	0.0					
	Actual pressure	Measured pressure	Indicated pressure		mmHg			1		
	60	59.90		240	1			[
	100	99.80			ł					
	140	139.80			Į					
	180	180.10		(BH 180	1					
	220	219.90		E	ł					
	260	259.80		Le (ţ		1			
	300			ssu	ł					
				JE 120	ţ					
			Run Cance							
				60	‡ /					
				0	¥					
					0	16	32	48	64	80
							Time (s)			

201.105.1 Long-term Automatic Mode

201.105 * Automatic cycling modes

201.105.1 LONG-TERM AUTOMATIC MODE

If an AUTOMATED SPHYGMOMANOMETER is equipped with a LONG-TERM AUTOMATIC MODE, a PROTECTION DEVICE shall be provided to ensure that:

- a) in NORMAL CONDITION:
 - the total duration of the alternating inflation/deflation periods in an unsuccessful DETERMINATION (see Figure 201.103) shall not exceed the maximum inflation time specified in 201.104; and
 - after each successful DETERMINATION;

the CUFF pressure shall be released and shall remain below the values in Table 201.102 for at least 30 s (see Figure 201.104); and

b) in SINGLE FAULT CONDITION:

if the duration of deflation below the values in Table 201.102 is less than 30 s (see Figure 201.105), then a pressure relief PROTECTION DEVICE functioning independently of the NORMAL CONDITION PROTECTION DEVICE, shall release the CUFF pressure to the values in Table 201.102.

Compliance is checked by functional testing.



Figure 201.104 - LONG-TERM AUTOMATIC MODE CUFF pressure in NORMAL CONDITION



2 Deflated time

CUFF pressure, Pc, as a function of time. NEONATAL MODE values in parentheses.

Blood Pressure Test

with Sys, Dias, HR and Two-end Sys/Dias Shift



3-2-6

Test Example

BPA700 Standard Assistant

201.105.1 Long-term automatic mode – Single Fault



3-2-6

Test Example

BPA700 Standard Assistant

201.105.1 Long-term automatic mode-Normal Condition



Advanced Specifications

Widen patient spectrum

Sphygmomanometer	Specs	BPA700
0-300	Manometer (mmHg)	0-400 ± (0.3% of reading + 0.5mmHg)
0-300	Pressure Source (mmHg)	20 to 400 ± 0.5mmHg
40-200(BPM), 40- <mark>290</mark> (PM)	Systolic range (mmHg)	25- <mark>300</mark>
40-195(BPM), <mark>30</mark> -245(PM)	Diastolic range (mmHg)	<mark>10</mark> -250
± 3	BP Accuracy (mmHg)	±1
40-200	Pulse rate (bpm)	30-300
±4% of reading	Pulse accuracy (bpm)	±1





Advanced Specifications

start of slow deflation

Widen patient spectrum

start of inflation



BP monitor



BPA700

Data Recording and Playback



- ▶ 3-2-8 Data Recording and Playback Process
- ▶ 3-2-9 Data Recording
- ▶ 3-2-10 Data Conversion
- ▶ 3-2-11 Data Playback

Blood Pressure Database Recorder and Player ISO 81060-2_5.2.3 REFERENCE DETERMINATION

3-2-8 Data Recording and Playback Process



Play : Data Playback and Algorithm Optimization









Record a Person's BP Raw Data

ISO 81060-2_5.2.4.3 Opposite arm simultaneous method

Note: Arteriosclerosis might affect measurements.



3-2-9 Data Recording

Record a Person's BP Raw Data

ISO 81060-2_5.2.4.1 Same arm simultaneous method





Convert BP Raw Data into Pulse Waveform



42

Play Oscillation Waveform of the Person

at Algorithm Development









Respiration Testing

- Respiration Module, paired with SECG 5.0 AIO, allows developers verifying BrPM in wider ranges and finer adjustments.
- Respiration Module, paired with SECG 5.0 AIO, makes possible RESP. Database Playback and cost-effectively enhances algorithms.



NIBP Testing

- P-T Graph of BPA700 GUI displays instant pressure changes over time by DUT.
- The Widest Ranges of SYS/DIAS(10-300mmHg) and Pressure Source (20-400mmHg).
- BPA700, also functioning as a BP Database Recorder and Player, <u>cost-effectively enhances</u> <u>algorithm testing</u> to ensure NIBP monitor cover a broader patient spectrum.
- Standard Assistant Software breaks down Test Requirements of Medical Standard into Pass Criteria, Preparation, Test Parameter per Condition, Auto-Verdict and P-T Graph, and also helps users keep records.





(4) Automated Testing Solution

Powered by WhaleTeq SDK



Key Features

✓ Free of Charge

Downloadable on the website, helping ramp up throughput for developers.

✓ User-Friendly

Well-documented APIs and Python sample codes available.

Cross-Platform Compatibility
 Supported OS: Windows, Linux and Raspberry Pi.

Example: IEC 60601-2-27 Test Sequence

Powered by SECG 5.0 AIO SDK

Test 201.12.1.101.11: Baseline Reset
print("Performing Baseline Reset test...")
secg.set_output_function(OutputFunction.Output_Sine)
secg.set_frequency(10) # 10Hz
secg.set_amplitude(1) # 1V
time.sleep(10)
print("Turning baseline reset off and measuring signal...")

Test 201.12.1.101.12: Pacemaker Pulse Display Capability
print("Performing Pacemaker Pulse Display Capability test...")
secg.set_output_function(OutputFunction.Output_Sine)
secg.set_frequency(1.5) # 1.5Hz
secg.set_amplitude(0) # 0mV
secg.set_pacing_rate(90) # 90bpm
secg.set_pacing_amplitude(700) # 700mV
secg.set_pacing_duration(2) # 2ms
time.sleep(10)
print("Measuring pacemaker pulse display...")

201.12.1.101.11 Baseline Reset

Signal Type, Frequency, and Amplitude Settings

201.12.1.101.12 Pacemaker Pulse Display Capability

Signal Type, Frequency, Amplitude, and Pacing Signal Settings

Example: IEC 80601-2-30 Test Sequence

Powered by BPA700 SDK

201.12.1.107

Reproducibility of Blood Pressure Determination

Test 201.12.1.107: Reproducibility of Blood Pressure Determination
print("Performing Reproducibility of Blood Pressure Determination test...")
for sample in ['A', 'B']:
 test_params.Systolic = 70
 test_params.Diastolic = 40
 test_params.HeartRate = 140
 for i in range(20):
 bpa.run_blood_pressure_test(test_params, None)
 time.sleep(10)
 print(f"Test {i + 1}/20 for sample {sample} completed.")

Systolic, Diastolic, and Heart Rate Settings for 20 Required Tests 201.12.1.107 Limits of the change in errorReproducibility of the BLOOD PRESSURE DETERMINATION

The laboratory-limits of the change in error reproducibility of the ELOOD PRESSURE DETERMINATION of the AUTOMATED SPHYGMOMANOMETER shall be less than or equal to 3.0 mmHg $(0.4\,\text{kPa})$.

Compliance is checked with the following test:

Two samples of the AUTOMATED SPHYGMOMANOMETER of the same MODEL or TYPE REFERENCE are needed to perform this test PROCEDURE.

NOTE. At the beginning of this compliance test network angule has been subjected to the mechanical stress tests of the general standard and the colleteral standards. Step Al subjects AUTOMATCD ENVIRONMANDEERER A to the stress tests and the taboratory limits of the change in error of the dLOOD PRESSURE DETERMINATION are compared before and attre these mechanical stresses.

a) Label one sample of the AUTOMATED SPHYGMOMANOMETER as A and the other sample as B.

a)D Prior to performing the other tests of this standard, adjust the PATIENT SIMULATOR to generate signals in such a way that the Autonates DEPARAMONETER displays approximately a pLASTOLIC BLOOD PRESSURE value of 40 mmHg (22-53 HPa) and a systolic slood PRESSURE value of 70 mmHg (48-003 JHPa) at a pulse rate of 140 beats/min in INEONATAL MODE and a DIASTOLIC BLOOD PRESSURE value of 80 mmHg (5-310.67 HPa) and a Systolic BLOOD PRESSURE value of 120 mmHg (30-716.0 HPa) at a pulse rate of 80 beats/min otherwise. Either sample of the AUTOMATED SEMPSGNOMMONDETER MAY be used for this step.

b)c) Perform 20 consecutive DETERMINATIONS with AUTOMATED SPHYOMOMANOMETER B. Calculate the means and standard deviations for both the DIASTOLIC BLOOD PRESSURE and the SYSTOLIC BLOOD PRESSURE.

c) Perform all the tests of this standard, except 201.106.

- d) Record these results as the AUTOMATED SPHYGMOMANOMETER B starting values.
- e) Verify that the standard deviation of the practicul succe pressure and of the systolic succe pressure are _2,0 mmHg (_s0,27 kPa) for the urtrearters privated and the systolic starting values [I either one of these criterion is not met, the combination of the simulator and AUTOMATED services and an another that the service of the
- f) Using the same PATIENT SIMULATOR and settings as in ab), perform <u>20 consecutive</u> DETERMINATIONS with AUTOMATED SPHYOMOMANOWETER A. Calculate the mean and standard deviation for both the pusatoric eloco PRESSURE and the systocic elocop PRESSURE.

g) Calculate the difference of the means calculated in b) and d) h) Ensure that the difference is below the limit.

- g) Record these results as the AUTOMATED SPHYGMOMANOMETER A starting values.
- h) Using AUTOMATED SPHYGMOMANOMETER Å, perform at least the following tests, without the simulation of single FAULT CONDITIONS, of this particular standard: 201.11.6.5, 201.12.1.102, 201.15.3.5.101, and 201.15.3.5.102 as well as IEC 60601-1:2005, 15.3.2, 15.3.3 and 15.3.4.
- Using the same PATIENT SIMULATOR and settings as in b), perform 20 DETERMINATIONS with AUTOMATED SHYVAMOMAKIONETER A. Calculate the means of the DIASTOLIC BLOOD PRESSUR and the SYSTOLIC BLOOD PRESSURE.
- j) Record these results as the AUTOMATED SPHYGMOMANOMETER A ending values.
- k) Using the same PATIENT SIMULATOR and settings as in b), perform 20 DETERMINATIONS with AUTOMATED SPHYSMOMANOMETER B. Calculate the means of the DIASTOLIC BLOOD PRESSURE and the systolic BLOOD PRESSURE.
- I) Record these results as the AUTOMATED SPHYGMOMANOMETER B ending values.
- m) For AUTOMATED SHYGMOMANOMETER B ending values, verify that the standard deviation of the DIASTOLIC BLOOD PRESSURE and of the SYSTOLIC BLOOD PRESSURE are ≤ 20 mmHg (< 0.27 KPG) if either one of these criterion is not met, the combination of the simulator and AUTOMATED SHYGMOMANOMETER has insufficient stability to perform this test PROCEDURE.
- n) For AUTONATED Service outware text B, verify that the absolute value of the difference between the mean starting values calculated in c) and ending values calculated in m) are = 2.0 mmHg (≤ 0.27 Ke), [if either one of these criterion is not met, the combination of the simulator and AUTONATED service wave the sinsufficient stability to perform this test processore.

b) For AUTOMATED SPHYGMOMANOMETER A, verify that the absolute value of the difference between the mean starting values calculated in f) and ending values calculated in i) are < 5,0 mmHg (< 0,07 KPa).</p>

① Prepare Sample A and Sample B of the same DUT model.
② B: 20* BP determinations → means & SD (B starting values)
③ A: 20* BP determinations → means & SD (A starting values)
④ Sample A alone goes thru mechanical stress test as per 201.12.1.102 & 201.15.3.5.101
(5) A: 20* BP determinations → means (A ending values)
⑥ B: 20* BP determinations → means & SD (B ending values)
⑦ Give verdicts as per criteria
 B ending standard deviation value: SD ≤ 2.0 mmHg (≤ 0.27 kPa)
• B abs. mean difference : ③ - ⑤ ≤ 2.0 mmHg (≤ 0.27 kPa)
 A abs. mean difference : 2 - 6 ≤ 5.0 mmHg (≤ 0.67 kPa)

Automated Testing Framework

Case Study

https://www.youtube.com/watch?v=LfxOEm_0n4o

ECG CMRR Test

ECG Performance Test

(5) Overall Conclusions

- Flexible ECG, PPG, Impedance RESP and NIBP Simulations At Once enhances analysis of vital sign correlation by patient monitors and verify exceptional spec. of developers.
- Replaying clinical databases and recorded raw data cost-effectively enhances algorithm testing.
- Standard Assistant Software breaks down and streamlines medical standard tests into Pass Criteria, Preparation, Test Parameter per Condition and Auto-verdict.
- SDK makes all WhaleTeq Products automatable for end-of-line testing.

Questions from Attendees

#	Questions
1	How can one evaluate the impact of changing the emitter brightness on the PPG signal quality?
2	Can BPA700 test smartwatches with inflatable strap?
3	Can PPG, coming from transflective sensor, also be tested?
4	Is ECG and IBP waveform synchronized? e.g. Change in ECG signal reflects change in IBP?

Q1: How can one evaluate the impact of changing the emitter brightness on

the PPG signal quality? TX Testing By SECG 5.0 AIO's PD

Equivalent to Sampling function on SECG SW

- Suggested steps - Use SECGEnableSampling with
SamplingCallback to extract
ADC counts of SECG 5.0 AIO's R&IR-PD.

2 Based on Developer's criteria, give verdict on DUT's LED performance.

Note: Developer may refer to the sample code in SDK available below.

Q2: Can BPA700 test smartwatches with inflatable strap?

Yes, but it would take one extra watch sample as an interface, which is illustrated as Watch A

59

Q2: Can BPA700 test smartwatches with inflatable strap?

Yes, but it would take one extra watch sample as an interface, which is illustrated as Watch A

a. Record

© 2025 WhaleTeq All Rights Reserve

b. Convert

© 2025 WhaleTeg All Risk Seserve

c. Replay

Q3: Can PPG, coming from transflective sensor, also be tested?

Yes, the PPG-2TF-660 can test it.

It's a transmittance PPG simulator, so it works for transmittance sensors.

PPG-2TF-660 (Transmittance SpO₂ Module) PPG-2TF-660 with the Multi Vital Sign Simulator SECG 5.0 AIO for ECG, PPG, Respiration, and PWTT Testing

Q4: Is ECG and IBP waveform synchronized? e.g. Change in ECG signal reflects change in IBP?

Actually IBP is not supported. We support NIBP testing.

Regarding IBP, we have a partner we've been working with for a long time and if you're interested, please let us know and we will direct you to the person.

Thank you for your time

Verify and Validate Your Product Design with Ease

service@whaleteq.com

+886-2-2517-6255 R)

WHALETEQ Facebook

