



【 From Theory to Practice 】
A Comprehensive Guide
to Conducting Accurate Tests for
ECG, PPG, PWTT, and Respiration Rate

Speaker: Richard Ding

Time	Contents
30 min.	<ul style="list-style-type: none">• ECG test methods according to international standards• Heart rate (HR) and SpO₂ testing using PPG technology• PWTT testing using ECG and PPG signals
20 min.	<ul style="list-style-type: none">• Respiration rate measurement using impedance and modulation test methods• Introduction of WhaleTeq's test equipment – SECG 5.0 AIO
10 min.	<ul style="list-style-type: none">• Conclusions• Q&A

ECG Test Methods

According to International Standards

ECG Standards

Three ECG Standards for Three Kinds of ECG

- **Diagnostic ECG: IEC60601-2-25:2011.**

A device intends to obtain, record and/or display a set of conventional or orthogonal ECG signatures, and provide a diagnostic ECG report. The tests include **Performance and database testing.**

- **ECG in patient monitor: IEC60601-2-27:2011.**

A monitor that obtains heart rate and waveform by ECG, is used to display the heart rate and/or heart rate waveform of the monitored patient, as well as an alarm for arrhythmia. The tests include **Performance testing only.**

- **Ambulatory ECG: IEC60601-2-47:2012.**

A system is capable of continuous recording and analysis or continuous analysis and partial or limited recording of an ECG, which records and stores the ECG and then analyzes it in a stand-alone unit, or records and analyzes simultaneously. The tests include **Performance and database testing.**

Performance Test Items

IEC60601-2-47

- Test ambulatory ECG **electrical performance**
- Five major performances for all ECG standards, i.e. **Amplitude**-related, **Input Impedance**, **Noise**-related, **Frequency Response**-related, and **Pacing Pulse**-related.
- Noise includes **mains frequency common mode noise** and **system (internal) noise**.
- Frequency response includes **low/high cutoff frequency response** and **pass band response**.

201.12.4.4.101	Linearity and dynamic range	←	Amplitude-related
201.12.4.4.102	Input impedance	←	Input Impedance
201.12.4.4.103	Common mode rejection	↙	Noise-related
201.12.4.4.104	GAIN accuracy		
201.12.4.4.105	GAIN stability	↘	Noise-related
201.12.4.4.106	System noise		
201.12.4.4.107	Multichannel crosstalk		
201.12.4.4.108	Frequency response	←	Frequency Response-related
201.12.4.4.109	Function in the presence of pacemaker pulses	←	Pacing Pulse-related
201.12.4.4.110	Timing accuracy		
201.12.4.4.111	GAIN settings and switching		
201.12.4.4.112	Temporal alignment		

Linearity and Dynamic Range – 201.12.4.4.101

Digital ECG

PURPOSE (DIGITAL ECG)


- Digital ambulatory recorders shall be capable of responding to and displaying an input signal of **10mV peak-to-valley (p-v)** in amplitude in the presence of a DC offset voltage of **$\pm 300\text{mV}$** .

METHOD

- Feed a **6.25Hz triangular wave, 0.5mV, 1mV, 2mV, and 10mV p-v** into the test circuit.
- Set switch S_3 to position **A** to superimpose an offset voltage of **$\pm 300\text{mV}$** .
- Alternatively:
 - **4Hz sine wave**, with the same amplitudes as above, either continuous or consisting of isolated cycles repeated once a second.

RESULTS

- ✓ Output signal amplitude referred to input shall not change by more than **10% or $50\mu\text{V}$** , whichever is greater.




Linearity and Dynamic Range Test

SECG 5.0 AIO Setup and Standard Assistant – Digital ECG

Signal Parameters		Signal add-on	
Signal Type	Triangle	DC Offset	ON <input checked="" type="checkbox"/> 300 mV
Frequency	6.25 Hz	Pacing	OFF <input type="checkbox"/>
	375 BPM	Noise	OFF <input type="checkbox"/>
Amplitude	10.00 mV	Sine wave add-on	OFF <input type="checkbox"/>
	100 ms	620kΩ / 4.7nF (on = short)	ON <input checked="" type="checkbox"/>
		Common mode to RL / N	OFF <input type="checkbox"/>

0.5mV, 1mV, 2mV, 10mV

±300 mV



$V_p-v \leq 10\%$ or $50\mu V$ whichever is greater

201.12.4.4.101 - Linearity and dynamic range Digital AMBULATORY RECORDERS			
Lead : Channel 1			
(mV)	DC Offset		
Amplitude	0 mV	300 mV	-300 mV
0.5 mV	0.50	0.50	0.50
1 mV	1.20	1.00	1.00
2 mV	2.00	2.00	2.00
10 mV	10.00	10.00	10.00

Input Impedance - 201.12.4.4.102

With a 620K Ω //4.7nF Parallel Circuit

PURPOSE


- The input impedance shall be **greater than 10M Ω** for the frequency specified in the test and for all input channels. This requirement shall be met across the total required **DC offset voltage of $\pm 300\text{mV}$** capabilities.

METHOD

- Apply a **10Hz sinusoidal, 5mV** amplitude.
- Connect the patient electrode connections of the first channel to P1 and P2. Connect all other patient electrode connections to P6.
- **Open S1 (620K Ω //4.7nF)** and measure the output amplitude change.
- Repeat the test with offset voltages of **300mV and -300mV** respectively.
- Repeat all these tests for **all other ECG channels**.


RESULTS

- The steady-state output amplitude shall not decrease by more than **6%**.



Input Impedance

Calculate by Voltage Divider Concept




$$V_2 = \frac{R_i}{R + R_i} * V_1$$


$$R_i = \frac{V_2}{V_1 - V_2} * R (620\text{K}\Omega)$$

If $V_2 = 0.8V_1$, $R_i = 4R \sim 2.5\text{M}\Omega$

If $V_2 = 0.94V_1$, $R_i \sim 16R \sim 10\text{ M}\Omega$



V1
= without $620\text{K}\Omega/4.7\text{nF}$



V2
= with $620\text{K}\Omega/4.7\text{nF}$
and $\pm 300\text{ mV DC}$


$V_2/V_1 \geq 94\%$

Input Impedance


SECG 5.0 Setup and Standard Assistant

Signal Parameters		Signal add-on	
Signal Type	Sine	DC Offset	ON <input checked="" type="checkbox"/> 300 mV
Frequency	10.00 Hz	Pacing	OFF <input type="checkbox"/>
	600 BPM	Baseline overload test	OFF <input type="checkbox"/>
Amplitude	5.00 mV	ECG special add-on	OFF <input type="checkbox"/>
	100 ms	620kΩ / 4.7nF (on = short)	OFF <input type="checkbox"/>
Pulse Width	100 ms	Common mode to RL / N	OFF <input type="checkbox"/>

S1 OFF: with 620KΩ
ON: without 620KΩ



V1
= without 620KΩ/4.7nF



$$V2/V1 \geq 94\%$$

V2
= with 620KΩ/4.7nF
and ±300 mV DC

201.12.4.4.102 - Input impedance	
(mV)	Patient Electrode
DC Offset	Channel 1
0mV(Vref)	4.90
0mV+620KΩ	4.70
300mV+620KΩ	4.60
-300mV+620KΩ	4.71
Report	Description Procedure Pass Criterion

Heart Rate (HR) and SpO₂ Testing Using PPG Technology

Green LED and PD to Measure HR


PPG: Photoplethysmography

LED emits green light to the skin


- ◆ **Partially absorbed** by skin tissue, arteries, and veins
- ◆ **Partially reflected** back to the PD of the sensor

Skin tissue and veins

- ◆ **In a static state**, the volume does not change with each heartbeat, and the intensity of absorbed light and reflected light is stable. It is converted into direct current (DC) by the PD.
- ◆ Darker skin color will absorb more incident light, so the reflected light will be weaker, and vice versa. So **the DC level represents the skin color range.**



- ◆ The converted electrical signal consists of two parts, **DC and AC.**
- ◆ The **AC signal comes with HR variation parameters.**




Artery

- ◆ **The blood volume changes** with each heartbeat and the heart rate, so it is converted into **alternating current (AC)** through PD.
- ◆ The more elastic the artery is, the larger amplitude changes the reflected AC signal has. So **the AC level represents the artery's elastic range.**

R & IR LED and PD

For R Value and SpO₂ Value of Pulse Oximetry

- Pulse oximetry determines SpO₂ value by illuminating vascular tissue with rapid switching between Red and IR light.
- AC signals of Red and IR PPG are sensitive to changes in SpO₂ value because of the variance in the light absorption of O₂Hb (Oxyhaemoglobin) and HHb (Deoxyhaemoglobin) at these two wavelengths.




- R value: using the amplitude ratio of AC/DC signals for both Red and IR wavelengths

$$R = \frac{(AC/DC)_R}{(AC/DC)_{IR}} \rightarrow SpO_2 = K1 + K2R \quad \text{R curve}$$

Note: SpO₂ value can be calculated as a linear function of R, where K1 and K2 are constants.


Heart Rate (HR) and SpO₂ Testing – Reflectance & Transmittance

Synchronization & 3 Major Parameters AC, DC, BPM, Effectively Simulating Reflection and Transmission Light




- **BPM** (Beat Per Minute) – AC signal's frequency for simulating different **heart rates**.
- **DC** – simulates brightness reflected of different **skin colors**.
- **AC** – simulates PPG waveforms changing in **arterial blood volume** within the skin.
- **PI (Perfusion Index)** – $(AC/DC) \times 100\%$

Reflectance



Transmittance




PWTT Testing Using ECG and PPG Signals

Pulse Wave Transit Time (PWTT)

Simultaneously Play ECG and PPG Waveforms

- Simultaneously play ECG and PPG waveforms, and **the time difference—PTTp (PTT peak) and PTTf (foot)**—of ECG and PPG waveforms is adjustable.
- The **estimated continuous blood pressure (BP) value** can be calculated by measuring the PTTp or PTTf.




Respiration Rate Measurement Using Impedance and Modulation Test Methods

Respiration Rate (RR) Measurement

Impedance Test Methods


- RR is defined as the measurement of a person's **breathing rate per minute (BrPM) at rest**.
- When the thoracic cavity expands/reduces (inspiration/expiration), **the chest impedance increases/decreases**.
- The ECG leads (Lead I, II, or III) can measure the transthoracic impedance changes through **a fixed current, high frequency signal source** generated from ECG internal respiration circuit.
- Test equipment uses **baseline and time-varying impedance circuits** to simulate the transthoracic impedance changes.
- The similar **baseline and time-varying impedance circuit (with different impedance ranges)** can be used to take the **ICG (Impedance Cardiography)** measurement.



Respiration Rate Measurement

Modulation Test Methods – EDR


- The method **extracting respiratory signals** via ECG or PPG technology is called **EDR (ECG-derived Respiration)**, and the other one using PPG technology is called **PDR (PPG-derived Respiration)**.
- **Three types of EDR** analysis to detect changes in respiration rate when breathing:
 - 1) The heart axis deviates and causes ECG amplitude to change (**AM: Amplitude Modulation**).
 - 2) The heart beating rate changes (**FM: Frequency Modulation**).
 - 3) The ribs and diaphragm position changes, which leads Electromyography (EMG) to change the ECG (**BM: Baseline Modulation**).



Respiration Rate Measurement

Modulation Test Methods – PDR

- **Three types of PDR** analysis to detect changes in respiration rate when breathing:
 - 1) The heart volume changes and causes the PPG amplitude to change (**AM: Amplitude Modulation**).
 - 2) The heart beating rate changes (**FM: Frequency Modulation**).
 - 3) The change in thoracic pressure causes changes in blood flow and affects the baseline change of PPG (**BM: Baseline Modulation**).



Introduction of WhaleTeq's Test Equipment – SECG 5.0 AIO

SECG 5.0 AIO

ECG & PPG Simulator for Performance Testing

- **ECG performance test** – built-in test circuits per defined in ECG standards and standard assistant
- **Heart Rate (HR) and SpO₂ testing** – Reflectance & Transmittance PPG modules
- **PWTT testing** – time difference between ECG and PPG signals
- **Respiration rate measurement** – Impedance Type, Wave Modulation Type (Baseline/Amplitude/Frequency)
- **Auto Sequence function** – assists in conducting preferred tests in sequence
- **Play raw data** – loading recorded or programmed waveform files facilitates verifying DUT algorithms
- **Software Development Kit (SDK)** – allows users to develop customized or automated test software



SECG 5.0 AIO

User Interface Introduction

ECG waveform

- Sine
- Triangle
- Square
- Rectangle pulse
- Triangle pulse
- Exponential
- ECG waveform**

Signal Parameters		Signal add-on	
Signal Type	ECG waveform	DC Offset	OFF <input type="checkbox"/> 0 mV <input type="checkbox"/> Variable
Frequency	1.00 Hz	Pacing	OFF <input type="checkbox"/>
		Noise	OFF <input type="checkbox"/>
		ECG special add-on <input checked="" type="checkbox"/>	
		620kΩ / 4.7nF (on = short)	ON <input checked="" type="checkbox"/>
Amplitude	2.00 mV	Common mode to RL / N	OFF <input type="checkbox"/>
Pulse Width	100 ms		

ECG special add-on parameters

AAMI EC 13 Drift test

Amplitude mV Frequency Hz

Respiration

Rate BrPM Apnea Time sec / min

Inhale:Exhale Ratio 1 :

Impedance

Basic Level Ω Variation Ω

Wave Modulation

Baseline % Amplitude % Frequency %

Time(s)

Special Functions

- IEC standard waveform
- File playback
- Add NST noise
- Auto Pacing / Heart Rate / Frequency Scan
- Lead Off
- Auto Sequence

Output Summary Dump Raw Data save SDK log save as script

Customized ECG Wave: 1.00 Hz (60 BPM), 2.00 mV

- PR interval 160 ms
- QT interval 350 ms
- P wave 0.20 mV
- P duration 100 ms

20 mm/mV 10 mm/mV 5 mm/mV 1 mm/mV 0.5 mm/mV 00:00:09

Amplitude (mV)

Time (s)


ECG **ECG Special Functions** Standard Assistant SpO₂ PWTT ▶ PLAY

Standard Assistant

- IEC 60601-2-25 / 27 / 47
- YY0782 / 0885 / 1079 / 1139
- GB9706.225 / 227
- YY9706.247

SECG 5.0 AIO

Respiration Rate Test



Respiration rate test spec.

BrPM Rate	0-170 BrPM in 1 BrPM step
Impedance Baseline	100 Ω ,200 Ω ,500 Ω -4500 Ω in 500 Ω step
Impedance Variation	0.0 Ω -5.0 Ω in 0.05 Ω step 5.0 Ω -10.0 in 0.1 Ω step

ICG test spec.

Impedance Baseline	20 Ω -100 Ω in 20 Ω step
Impedance Variation	0.0 Ω -1.0 Ω in 0.05 Ω step 1.0 Ω -5.0 in 0.1 Ω step

SECG 5.0 AIO

PWTT Test

PWTT

ECG

Load ECG

BPM

Amplitude mV

T Amplitude mV

QRS Duration ms

SpO₂

R IR Load from SpO₂ page

SpO₂ % Sync Pulse

PI % Inverted

DC mV Lock DC

AC mV Lock AC


Time Difference

PTTp ms Lock PTTp

PTTf ms Lock PTTf

Ambient Light

Respiration



The image shows a screenshot of the SECG 5.0 AIO PWTT Test interface. On the left, there is a graph displaying an ECG waveform (blue line) and a corresponding SpO₂ waveform (red and purple lines). A horizontal double-headed arrow labeled 'PTTp' indicates the time difference between the R-peak of the ECG and the start of the SpO₂ pulse. The interface is divided into several control panels: 'ECG' with settings for Load ECG, BPM (60), Amplitude (1.0 mV), T Amplitude (0.20 mV), and QRS Duration (100 ms); 'SpO₂' with settings for R (checked), IR, Load from SpO₂ page, SpO₂ (100%), PI (2.000%), DC (625 mV), AC (12.50 mV), Sync Pulse (checked), Inverted, Lock DC (checked), and Lock AC; 'Time Difference' with settings for PTTp (350 ms) and PTTf (200 ms), both with lock options; and 'Ambient Light' and 'Respiration' both set to 'OFF'.

© 2022 WhaleTeq All Rights Reserved

25

Conclusions

ECG, PPG, PWTT, RR, ICG and Test Equipment

- There are three ECG standards for three kinds of ECG performance and database tests.
- Five major performances for all ECG standards, i.e. **Amplitude**-related, **Input Impedance**, **Noise**-related, **Frequency Response**-related, and **Pacing Pulse**-related.
- PPG technology is currently applied mainly to Heart Rate (HR) and SpO₂ measurements.
- Three major parameters: **BPM** as HR, **AC** level represents the artery's elastic range, and **DC** level represents the skin color range.
- The **PWTT** can be adjusted to effectively verify the accuracy of the blood pressure measurement algorithm.
- **RR and ICG** (Impedance Cardiography) measure skin impedance changes to get **BrPM and CO** (Cardiac Output).
- Test equipment simulates the impedance changes with **baseline and time-varying impedance circuits**.
- **SECG5.0 AIO** can support **ECG, PPG, PWTT, RR, and ICG** applications.



Verify and Validate Your Product Design with Ease



www.whaleteq.com



service@whaleteq.com



+886-2-2517-6255



WHALETEQ
Web



WHALETEQ
Facebook



WHALETEQ
YouTube



WHALETEQ
LinkedIn