

GS Elektromed. Geräte G. Stemple GmbH · Hauswiesenstr. 26 · 86916 Kaufering / Germany

WHALETEQ CO., LTD
9F-1, NO. 104, Minquan W. Rd.,
Taipei City 10361, Taiwan

**GS Elektromedizinische Geräte
G. Stemple GmbH**
Hauswiesenstraße 26
86916 Kaufering / Germany
Tel. +49 8191 65722-0
Fax +49 8191 65722-22
info@corpuls.com
www.corpuls.com

Ihr Zeichen · Your reference

Unser Zeichen · Our reference

Direktwahl · Direct dialing

E-Mail · e-mail

Datum · Date
01.04.2017

Subject: Letter of Recommendation for WhaleTeq SECG SDK (Software Development Kit)

Dear WhaleTeq team,

I am writing you this letter as a letter of recommendation to the WhaleTeq SECG SDK provided by your company.

Your WhaleTeq SECG SDK enables us to integrate the WhaleTeq SECG Test Unit Set into our automated hardware-in-the-loop test setup. With this test setup we can not only send ECG signals from a PC to our product but rather get the response of our product back to our PC and compare them with an expected result. This way we can automatically perform a variety of test cases after each adaption of our algorithms and validate our products.

It is important to validate the overall buildup in order to use this test setup for our mandatory standard test. Concerning this matter the validation of the WhaletEQ SECG Test Unit Set towards the standards 60601-2-25:2011, 60601-2-27:2011, 60601-247:2012, 60601-2-47:2012 and ANSI/AAMI EC 13 saves us a lot of time and efforts.

Therefore I can recommend the WhaleTeq SECG SDK as a helpful tool for the validation of ECG products.

Sincerely,



Dr. Antoun Khawaja

Director signal processing

GS Elektromedizinische Geräte G. Stemple GmbH