



WhaleTeq CMRR series PC Software

Version: 1.0.7.3 Build: 20191218

=====

1. for CMRR3.0E , remove all assistant tests except for IEC 80601-2-26

Version: 1.0.6.9 Build: 20190828

=====

1. fix incorrect DC offset text in the Monitor
2. CMRR3.0E : update IEC-2-26 assistant name
3. fix typo

Version: 1.0.5.11 Build: 20190715

=====

1. added support DC Offset +/- 300 to V2~V6 in all standards (except NOISE)
2. set a default setting when device is connected (Standard: "Manual", Voltage: "Off", Frequency: "50 Hz", Electrode: "None with impedance", DC Offset: "Off")

Version: 1.0.4.2 Build: 20180907

=====

1. Update ui layout
2. Update IEC 2-26 test description

Version: 1.0.3.2 Build: 20171003

=====

1. Added support to check the latest AP version and upgrade
2. Added support to update the activation key (in the system menu -> 'License...')

Note:

1. CMRR series including CMRR3.0 & CMRR3.0E.