

**One-Stop Solution**



# Comprehensive IEC 80601-2-30 Testing Package

World's first ever complete test solution for NIBP blood pressure monitor with waveform record and playback function, simplifying and reducing the cost of clinical trials.

## Simulator & Software

- ① **BPA700** Blood Pressure Simulator & Analyzer
- ② **Envelope Adjustment Software**
- ③ **IEC80601-2-30 Standard Assistant Software**

### Accessories

- ④ T-connector x 1
- ⑤ Connecting Tubes (Port Side) x 5
- ⑥ 100ml Air Reservoir x 1
- ⑦ 500ml Air Reservoir x 1

### Training & Service

- BPA700 Operation Training (2 hours)
- IEC 80601-2-30 Testing Training (2 hours)
- Testing Consultation (1 hours)

**Optional Items** • Waveform Playback Software • C3R3 3-Year Calibration Service and Warranty Extension

Package Contents

# BPA700 Installation – NIBP Monitor Verification

Simulate various blood pressure conditions observed in real clinical settings.

Able to connect a blood pressure monitors directly to patient to record and store patients' pressure values for data playback, facilitating validation and optimization of the monitor's algorithms.

