



**Test Solutions for**  
Medical Device Manufacturers

# AIMS

## Bring Your AED Readiness and Compliance

A complete AED after-sales solution to manage your AED programs with peace of mind. Get your AEDs ready in three steps: regular maintenance test, real-time monitoring, and data management, helping you lower the legal risk of AED ownership.

- First AED after-sales maintenance system designed for service partner with simple, affordable and complete AED testing solutions
- Unified AED Management: Up-to-date reminder, test report, real-time status tracking and complete data management
- Complete Testing Solution: Inspection, shock energy test, and independent battery test with compliance
- Multi-Device Monitoring: Get push notifications by message or email in a timely manner



**AIMS**  
Scan for  
more information

 [www.whaleteq.com](http://www.whaleteq.com)

 [service@whaleteq.com](mailto:service@whaleteq.com)

 8F., No. 125 Songjiang Rd., Zhongshan Dist., Taipei City 104474, Taiwan

 +886-2-2517-6255

 +886-2-2596-0702

## • Specifications

Module	Description
AED CMMS PLATFORM	Monitoring Testing Information <ul style="list-style-type: none"> <li>• Use the indicator light to display the alarm status of the device, including operating status, device data and monitoring records</li> <li>• Display on regional map or as a table with filtering function</li> </ul>
	Task Management <ul style="list-style-type: none"> <li>• Five major types of work orders: Installation, alarm, maintenance, training, and other orders</li> <li>• Classified content and status of each work order</li> <li>• Email notification of real-time work order status</li> <li>• Easy-to-use filter with user-defined criteria</li> </ul>
	Sales Management <ul style="list-style-type: none"> <li>• Common brands AED models &amp; accessories are built-in as model list</li> <li>• Adopt progressive hierarchical management for AED, including inventory information and sales records</li> </ul>
	Report Management <ul style="list-style-type: none"> <li>• The interface has modification, submission, sign-off, and transmission, and viewed inspection details</li> <li>• Report supports multi-language and logo substitution</li> <li>• Report summaries health status of AED and supplements IEC 60601-2-4 recommendations</li> </ul>
	Statistical Analysis <ul style="list-style-type: none"> <li>• Statistical Data for Sales, work order and the AED monitoring status</li> </ul>
DFS200 APP	Equipment Management <ul style="list-style-type: none"> <li>• List complete equipment, display equipment status through different indicators</li> </ul>
	Order Management <ul style="list-style-type: none"> <li>• Arrange the route of the work order to be processed</li> <li>• View the content and details of the work order</li> <li>• Upload on-site photo records</li> <li>• After the work order is confirmed, sign and archive records</li> <li>• Work Order Report with complete test results</li> </ul>
	Manual Testing <ul style="list-style-type: none"> <li>• Set the test signal output and pass-criteria</li> <li>• Various test sequence (up to six tested waveforms) and test file corresponding to different device brands / models are available for customization</li> <li>• Display the test results and output waveform from AED</li> </ul>
	Battery Testing <ul style="list-style-type: none"> <li>• Use probe with load to measure battery voltage and estimate remaining battery level</li> <li>• Data record as long-term trend observation</li> </ul>

## • Ordering Information

### ♦ Test System

Module	Description
DFS 200 AED Field Tester Set	The set contains DFS200 AED Field Tester, its standard accessories, USB Cable, Carrying Bag, and Battery Measurement Probe with Load
License Fee for AED CMMS PLATFORM_MOQ 100 up	AED CMMS PLATFORM Channel License Fee for MOQ100 and above
License Fee for AED CMMS PLATFORM_MOQ below 100	AED CMMS PLATFORM Channel License Fee for MOQ less than 100 (Minimum order 10 and above)

### ♦ Test System Optional

Module	Description
Service Charge Rental for DFS200	Monthly Rental Fee of DFS200 AED Field Tester, its standard accessories, Carrying Bag, and Battery Measurement Probe with Load

## • Test System Structure

1. Testing APP    2. DFS200    3. AED CMMS    4. IoT Module (Customization)

