

# WhaleTeq DFS400 SDK Changelog

---

## Version 1.0.2.0 - 2024-07-15

---

1. NEW: support Linux platform

## Version 1.0.1.1 - 2024-05-09

---

1. FIX: read/write license on wrong address for DFS360

## Version 1.0.1.0 - 2024-05-09

---

1. CHG: merge DFS360 in to DFS400 SDK
2. NEW: support play data

## Version 1.0.0.14 - 2023-09-18

---

1. CHG: NSR BPM min/max

## Version 1.0.0.12 - 2023-06-02

---

1. CHG: get pulse buffer size by command

### Note:

1. only compatible with FW version 0.42.X

## Version 1.0.0.11 - 2023-05-26

---

1. CHG: get/clear energy computed flag by command

## Version 1.0.0.10 - 2023-03-30

---

1. CHG: can't set load immediately after receive shock
2. CHG: can't set any parameters when ECG is playing
3. CHG: NSR sample rate 1000 to 500

## Version 1.0.0.9 - 2023-03-29

---

1. NEW: set load will sleep 1 second when switching under 1 seconds

**Note:**

1. only compatible with FW version 0.40.1(DFS400)/0.40(DFS360) and later

## **Version 1.0.0.5 - 2022-12-14**

---

1. CHG: connect device in demo mode can play ECG

## **Version 1.0.0.4 - 2022-12-13**

---

1. FIX: connect check demo mode

## **Version 1.0.0.3 - 2022-12-08**

---

1. FIX: NSR frequency and amplitude do not update
2. NEW: can connect DFS 360/400
3. NEW: energy ratio change to charge ratio
4. NEW: set load resistance and play ECG sample (Python)

## **Version 1.0.0.2 - 2022-11-02**

---

1. NEW: Set DAC value command
2. NEW: Support demo mode
3. NEW: Get serial number

## **Version 1.0.0.1 - 2022-10-17**

---

1. initial