

# WhaleTeq WECG400 Changelog

---

## Version 1.0.2.1 - 2025-07-14

---

1. CHG: set DC Offset to 0mV after closing the application

## Version 1.0.1.8 - 2024-04-26

---

1. CHG: disable functions which not available in previous generation
2. FIX: after the device is connected, the input impedance is not selected to a default item

## Version 1.0.1.6 - 2024-03-08

---

1. FIX: failed to set input impedance 200 M $\Omega$  / 10 pF

## Version 1.0.1.5 - 2024-02-22

---

1. NEW: added support to configure the connection between Signal Ground and GND Ground; only compatible with fw 1.4.32 and later
2. FIX: input impedance - 10 M $\Omega$  / 180 pF
3. NEW: added input impedance - 200 M $\Omega$  / 10 pF (requires hardware support)

## Version 1.0.1.1 - 2023-05-10

---

1. CHG: enhance error handling for ESD

## Version 1.0.0.8 - 2022-09-23

---

1. CHG: update command protocol; only compatible with fw 1.4.25 and later

## Version 1.0.0.7 - 2022-09-06

---

1. FIX: blocking issue if the length of the raw data is too short
2. FIX: failed to run the minimum frequency of ECG to 0.02Hz
3. FIX: incorrectly play raw data

## **Version 1.0.0.3 - 2022-09-02**

---

1. FIX: not accurate rectangle pulse width
2. FIX: incorrect pulse width range

## **Version 1.0.0.2 - 2022-08-25**

---

1. FIX: failed to set amplitude

## **Version 1.0.0.1 - 2022-08-18**

---

1. initial