

WhaleTeq SEEG 100/100E Changelog

Version 1.0.10.1 - 2024-06-01

1. CHG: update the download link of eegmmidb database
2. CHG: .NET Framework 4.5.2 is required

Version 1.0.9.1 - 2023-12-22

1. FIX: when playing a file, a repeated waveform is generated in the beginning

Version 1.0.8.5 - 2023-12-04

1. FIX: assistant name typo
2. FIX: the pulse waveforms cannot reach 5Hz
3. CHG: remove date information when playing special waveforms
4. FIX: in simulation mode, the time keeps updating after stopping special waveform output
5. FIX: when outputting the signal, the signal will continue to be delayed

Version 1.0.7.1 - 2023-06-26

1. FIX: failed to run if the device is not ISP version

Version 1.0.6.6 - 2023-06-21

1. FIX: incorrect check version of live update

Version 1.0.6.5 - 2023-06-07

1. FIX: incorrectly set the output lead when running assistant

Version 1.0.6.2 - 2022-12-12

1. CHG: support assistant GB9706-226.

Version 1.0.6.1 - 2021-09-22

1. FIX: the output amplitude in the graph of the main window is incorrect when SEEG 100E is connected

Version 1.0.5.10 - 2021-05-19

1. Fix "Under non 100% Windows resolution, main form would shrink after assistant form shown" problem.

Version 1.0.5.9 - 2019-08-02

1. Support SEEG-100E(Assistant 80601-2-26:2019).

Version 1.0.5.8 - 2019-03-13

1. Fix first time database downloading failed bug.
2. Change limit of "Variable DC Offset" in Calibration from 5% to 1%.
3. Fix "Live Firmware Update" notification failed bug.
4. Fix Serial Number display error on windows bar.

Version 1.0.5.6 - 2018-03-09

1. Live Update for Software and Firmware.
2. Change the initial value of amplitude and frequency to "100uV" and "10Hz".
3. Refine "Load EDF Database" and "Raw Data Replay".
4. Installation version with Certification supported.
5. Support "Epileptic Database" -
chbmit(<https://physionet.org/physiobank/database/chbmit/>)

Version 1.0.0.3 - 2017-02-09

Add ISP (In System Programming) Feature, it includes below two functions:

- Firmware Update
- Binding Assistant active key to SEEG device. The feature is only supported with new hardware and firmware.

Version 1.0.0.2 - 2016-05-19

1. New Feature: Support EDF file and File Manager

Version 1.0.0.1 - 2015-10-28

1. Initial release version for EEG AP.