

WhaleTeq CMRR 3.0 Changelog

Version 1.0.10.2 - 2024-07-16

1. NEW: added support FW version 4.x

Version 1.0.9.7 - 2024-04-16

1. NEW: added support GB 9706.225/ GB 9706.227/ YY 9706.247/ GB9706.226

Version 1.0.8.1 - 2022-07-13

1. FIX: (CMRR3.0E) set the impedance option from "All" to "None" if standard NOISE is selected

Version 1.0.7.3 - 2019-12-18

1. CHG: (CMRR3.0E) remove all assistant tests except for IEC 80601-2-26

Version 1.0.6.9 - 2019-08-28

1. FIX: incorrect DC offset text in the Monitor
2. CHG: (CMRR3.0E) update IEC-2-26 assistant name
3. FIX: fix typo

Version 1.0.5.11 - 2019-07-15

1. CHG: added support DC Offset +/- 300 to V2~V6 in all standards (except NOISE)
2. CGH: set a default setting when device is connected (Standard: "Manual", Voltage: "Off", Frequency: "50 Hz", Electrode: "None with impedance", DC Offset: "Off")

Version 1.0.4.2 - 2018-09-07

1. CHG: update ui layout
2. CHG: update IEC 2-26 test description

Version 1.0.3.2 - 2017-10-03

1. CHG: added support to check the latest AP version and upgrade
2. CHG: added support to update the activation key (in the system menu -> 'License...')